


Thermo Scientific Microm HM 560 Cryostat-Series



Designed to be above

The perfect solution for your cryo application



*Accessories not included. **Cryostat stage optional.

Unachieved performance for all applications in cryotomy

The cryostat HM 560 sets new standards with its combination of highest operator comfort and reliability with reproducible, optimal sectioning results. With this new cryostat line, Microm achieved to combine ergonomics and user protection with the renowned efficiency of our instruments. The well-known modularity with the options Vacutome, specimen cooling, cutting drive motor, active deep freezing device as well as the height adjustment offers for each application and individual need the optimal instrument concept.

Unique freezing system by controlled specimen and knife cooling

Due to the controlled knife temperature, the later rolling up of the section below the anti-roll guide is avoided. The normal heating up of the blade and/or knife carrier during ribboning with an open sliding window can be avoided as well. As the temperature of specimen and knife and/or blade carrier can be set independently of each other, you have the opportunity to adjust them to your tissue in an optimal way. The specimen cooling is integrated into the microtome. This way, the chamber remains in its entire volume.

Efficient fast freezing station to avoid freezing artefacts

An integrated Peltier element allows a surface temperature of down to -60°C resulting in a fast freezing-on of the tissue onto the chuck. Freezing artefacts can mainly be avoided via these deep temperatures.

Tissue detach function inside the cryochamber

The temperature of $+5^{\circ}\text{C}$ on the integrated and temperature controlled detach element guarantees a fast detaching inside the cryochamber which helps avoiding carryover of pathogens from the cryochamber.

Precision and constancy of the section thicknesses

A precise stepping motor for the feed of the selected section thickness guarantees constant and reproducible section thicknesses. The section thickness is set by pressing a key outside the chamber. Trim and fine section thickness are both shown on the display. The section thicknesses can be set from $0,5\ \mu\text{m}$ to $500\ \mu\text{m}$.

Highest operator comfort and maximum user benefit

The ideal positioning of all operating elements allows for an easy operation. The desired values of knife and/or specimen temperatures, trim and fine section thicknesses, cutting window, activation of the automatic approach and fast freezing station can easily be activated via the keyboard. Both the set values and the actual temperatures are shown on the large LC display. Ideal temperature settings for different tissue samples can be stored. Adjustment of cutting window for disposal and stretching of the section via the Vacutome as well as for the motorized cutting drive is carried out by pressing a button.



Electronical handwheel brake

The electronical handwheel brake can be adjusted in any position.

Innovative knife carrier for conventional knives and disposable blades

This patented knife carrier can be used for both conventional knives and disposable blades by a simple alteration. The anti-roll guide can be oriented optimally to the cutting edge of the blade and/or knife.

Adjustable chamber illumination and heated arm rests

The chamber illumination which is integrated in the handle of the sliding window can be set to the desired position to avoid reflections onto the anti-roll plate and to illuminate the cryochamber optimally. A heated arm rest increases the comfort.

Precise specimen orientation

The specimen can be oriented 8° in either direction via the specimen fine orientation. The parallel positioning of the chuck towards the knife is guaranteed via the zero position. The "quick release" mechanism for clamping chucks allows for a fast changing of specimens.

Protection against debris inside the cryochamber

The spacious cryochamber is made of stainless steel which can easily be cleaned. The rotary microtome is completely outside the chamber and thus protected against debris. The easily removable section waste tray with integrated brush shelf is positioned below the blade carrier and/or specimen.

Modularity for individual requirements

Vacutome: Section stretching and decontamination of section waste

The Vacutome is an integrated, patented module for stretching sections already during sectioning. Thin sections can be generated more easily. Section waste is automatically transferred into a coarse filter. The suction air is then filtered via a micro filter.

Motorized cutting drive

The motorized cutting drive is a valuable feature for routine work as well as for problematic tissue. The cutting speed can be selected continuously depending on the respective needs.

Ergonomy for sophisticated needs

The working height of the cryostat can be adjusted motorized and continuously to your desired working position or size. A storage surface on the upper side of the cryostat hood offers enough space for staining kits.



Technical data/accessories



Cryostat stage with height adjustment

Depending on the version, the cryostat can be adjusted to the desired working position either mechanically or motor-driven. The vertical adjustment range is 180 mm.

Neg -50 Cryomedium

Using the freezing medium Neg -50, the automatic approach can be carried out down to a specimen temperature of -50°C.

Disposable blades

SEC 35 blades are made of high-quality steel and coated with a unique, newly developed layer allowing for longevity of the knife at best cutting results.

Additional storage unit for chucks

Inside the cryochamber, this accessory serves as an additional storage unit for six chucks. Outside the cryochamber, the storage unit is used for detaching the tissue specimens from the chuck.

Technical data HM 560

Indication of data such as set and actual temperature, section sum, remaining travel, number of sections, actual time, trim and fine sectioning thicknesses, retraction function

Automatic approach system (ACA)

Electronical handwheel brake – can be locked in any position

Separate temperature control for knife carrier down to -35°C

Temperature for specimen clamping +10°C to -50°C

Heated sliding window

Standby and sleep mode for saving energy

Automatic defrosting device with immediate defrosting and additional interruption function

Detach function for the fast detaching of specimens from the chucks

Heavy-duty rotary microtome

Section thicknesses from 0,5 µm to 100 µm:	Trimming section thicknesses from 5 µm to 500 µm:
up to 2 µm in 0,5 µm-increments	5 – 10 µm in 5 µm-increments
up to 10 µm in 1 µm-increments	10 – 100 µm in 10 µm-increments
up to 20 µm in 2 µm-increments	100 – 200 µm in 20 µm-increments
up to 50 µm in 5 µm-increments	200 – 500 µm in 50 µm-increments
up to 100 µm in 10 µm-increments	

Specimen retraction that can be turned off

Horizontal feed range: 48 mm

End position indication and automatic turn-off when reaching the front and rear end positions of the horizontal feed.

Motorized coarse feed and reset movement with three different speed settings

Vertical cutting stroke: 60 mm

Max. specimen size: 70 x 55 mm

Specimen can be rotated in Z-axis, 360° and specimen orientation

0-positioning of the chuck

Cutting speed can be controlled continuously

Accelerated return travel for fast operation

Cutting window that can be selected

Optional

Three operating modes: interval, single and continuous stroke

Can be activated either via the foot pedal or manually

Motor brake and emergency stop switch for safe operation

Integrated Vacutome system for stretching cryosections as well as for the suction of debris

Integrated fast freezing device with separate cooling element (Peltier) down to max. -60°C

Dimensions: (W x D x H): 769 (w/o handwheel) x 932 x 578 (w/o stage) mm
Weight: 141-145 kg (depending on the model)



Rev. 11/08

© 2009 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Anatomical Pathology

4481 Campus Drive
Kalamazoo, MI 49008
(800) 522-7270

Robert-Bosch-Str. 49
69190 Walldorf/Germany
+49(0)6227-8360

93/96 Chandwick Road
Astmoor, Runcorn Cheshire WA7 1 PR
0800 018 9396, + 1 (0) 1928 562501

www.thermo.com/pathology

Thermo
SCIENTIFIC