




Incubation box

SlideFix



The SlideFix is a range of 3 purpose built humid incubation boxes with special polymer rails that grip your slides, preventing sample loss.

It has proved to be the ideal system to semi-automate immunohistochemistry methods, we are confident that once you use it, you will wonder how you managed before.

clear or black lids

tap or hose can be added

stackable

range of 3 sizes

range -10°C to +45°C

ideal as a slide workstation

8 years of sale

All slide show base units are made from ABS, lids from ABS or PETG and are able to withstand a wide variety of chemicals, but not Chlorinated Hydrocarbons.

All base units and lids are made from a single sheet which is heat formed into shape, this provides greater strength.

SlideFix 10 – Holds up to 10 slides

Unit with clear lid - cat no 582 430 510
Unit with black lid - cat no 582 430 511
Clear lid only - cat no 582 430 512
Black lid only - cat no 582 430 513

SlideFix 20 – Holds up to 20 slides

Unit with clear lid - cat no 582 430 520
Unit with black lid - cat no 582 430 521
Clear lid only - cat no 582 430 522
Black lid only - cat no 582 430 523

SlideFix 30 – Holds up to 30 slides

Unit with clear lid - cat no 582 430 530
Unit with black lid - cat no 582 430 531
Clear lid only - cat no 582 430 532
Black lid only - cat no 582 430 533

Cellab Nordia AB
08-344 770
www.cellab.se