



## See how easy it can be with Thermo Scientific EZ Cytofunnels®

The Thermo Scientific Cytospin technology has been synonymous with cyto centrifugation for over 40 years. It provides economical thin-layer preparations from any liquid sample. Using centrifugal force, the cyto centrifuge deposits cells onto a clearly defined area of a glass slide while residual fluid is absorbed into the cytofunnel's filter card. Choose EZ funnels for all your cyto centrifugations needs and simplify your work.

The Cytospin 4 and EZ Cytofunnels are the most cost-effective way to produce thin-layer preparations. Safe, easy to use, and reliable, this method optimizes cell yield while improving laboratory efficiency.

Learn more at [www.thermo.com/pathology](http://www.thermo.com/pathology)  
800-522-7270

### No More Metal Clips!



EZ funnels eliminate the need for metal clips and are disposable. The EZ funnels improve turn-around time and produce the most cost-effective, thin-layer preparations available today.



See the reverse side for a complete list of available products and ordering information

# Thermo Scientific EZ Cytofunnel Sample Chambers

For convenient disposable  
economical use.



## EZ Cytofunnel Features:

- Single use for all Thermo Scientific Cytospin® Cytocentrifuge applications
- Optimizes cell yield
- Reduces the risk of exposure to pathogenic samples
- Disposable to eliminate time-consuming decontamination
- Decreases the possibility of cross-contamination
- Eliminates the need for metal Thermo Scientific Cytoclips®
- Produces more cost-effective, thin-layer preparations when compared to other thin-layer methods.

## Includes:

Closure caps and easy-to-use operational illustrations



## EZ Single Cytofunnel with White Filter Card

- For use with sample volumes up to 0.5ml
- Cell deposition area of 6mm in diameter (28mm<sup>2</sup>)
- Pre-attached white filter card for high absorbency

Qty.	Ordering Number
40/box	A78710003
200/pk	A78710020



## EZ Single Cytofunnel with Brown Filter Card

- For use with low sample volumes of 0.4ml or less
- Cell deposition area 6mm in diameter (28mm<sup>2</sup>)
- Pre-attached brown filter card for slower absorbency
- Excellent for hypocellular specimens such as CSF

Qty.	Ordering Number
40/box	A78710004



## EZ Single Cytofunnel Bench Stand

- Convenient for loading specimen samples outside of the Cytospin chamber
- Stainless-steel construction, autoclavable
- Holds 12 EZ Cytofunnels

Qty.	Ordering Number
Ea.	A78710022



## EZ Double Cytofunnel

- Accommodates two sample deposition areas on a single slide
- Specimen samples are easily identified on pre-marked slides
- Pre-attached filter card
- Ideal for use with immunochemistry techniques

Qty.	Ordering Number
40/box	A78710005
200/cs	A78710021



## EZ Megafunnel™ with Cap and Slide

- For use with samples up to 6ml
- Minimizes the need to pre-centrifuge
- Cell deposition area of 22 x 14.75mm (325mm<sup>2</sup>) for increased cell yield
- Suitable for both air-dried and fixed preparations, providing additional valuable diagnostic information

## Includes

Thermo Scientific Megafunnel Cytoslides®

Qty.	Ordering Number
25/box	A78710001