


Thermo Scientific SlideMate On-Demand Slide Printer

Reduce Specimen Identification Errors in Cytology Slide Processing

Reduce the risk of loading the incorrect slide/vial combination into the Cytology slide processor by printing slides on-demand with the SlideMate slide printer. The small size of the SlideMate slide printer brings the convenience of on-demand slide printing directly to the slide processor, allowing the user to print one slide at a time for a single vial. The result of this process is that only one slide— the correct one— is available to load with the vial.

- **Small footprint fits easily into the Cytology preparation area**
- **On-demand slide printing reduces errors and improves productivity**
- **Prints directly onto the slide eliminating the need for labels or handwriting**
- **Chemical resistant media ensures print reliability**



“The elimination of manually applying labels and hand writing patient names on slides by using the SlideMate slide printer has virtually eliminated human error with this critical step in the cytology process. Additionally, the SlideMate slide printer has paid for itself in 12 months since we started using it last year due to discontinuing purchases of expensive xylene resistant slide labels.”

Francesco Mercuri, MS CT(ASCP), MT(ASCP)
Cytology/Histology Manager – Enzo Clinical Labs

“Printing slides using SlideMate™ slide printer has reduced potentially serious patient safety errors, allowed us to lean our process and improve workflow.”

Paul Miller
Anatomical Pathology Manager,
Children’s Hospital of Philadelphia



Ordering Information

Description	Order Number
SlideMate Print On-demand Printing system <i>Includes SlideMate Printer, scanner, keyboard interface, keyboard</i>	B81300001
SlideMate with LabWriter software <i>Includes SlideMate Printer, LabWriter Software for SlideMate</i>	B81200002
SlideMate with PC – Premium Package <i>Includes SlideMate Printer, All-in-one PC, digital imager, LabWriter software</i>	B81300003
SlideMate Printer	B81300004
SlideMate with Mini PC (built-in digital imager, clamp style arm) – Premium Package <i>Includes SlideMate printer, mini all-in-one PC with built-in digital imager and clamp style arm, LabWriter Software, cross-over cable, keyboard, mouse, USB hub, stylus</i>	B81300016
SlideMate with Mini PC (built-in digital imager, desk-top stand) – Premium Package <i>Includes SlideMate printer, mini all-in-one PC with built-in digital imager and desktop stand, LabWriter Software, cross-over cable, keyboard, mouse, USB hub, stylus</i>	B81300017
SlideMate with Mini PC (clamp style arm) – Premium Package <i>Includes SlideMate printer, mini all-in-one PC with clamp style arm, LabWriter Software, cross-over cable, hand-held digital imager with stand, keyboard, mouse, USB hub, stylus</i>	B81300018
SlideMate with Mini PC (desktop stand) – Premium Package <i>Includes SlideMate printer, mini all-in-one PC with desktop stand, LabWriter Software, cross-over cable, hand-held digital imager with stand, keyboard, mouse, USB hub, stylus</i>	B81300019

Specifications

Width	6.0 in, 152mm
Height	8.9 in, 226mm
Depth	10.4 in, 264mm
Weight:	9.5lbs., 4.3 kg.,
Power Requirements:	Universal 100 – 240 VAC, 47 – 63 Hz, 1.2A Max.
Color	Sterile Steel
Slide Type	3in x 1in or 25mm x 75mm x 1mm with Superfrost or Colorfrost end
Slide Capacity	150 pieces (input stack) 10 pieces (output tray)
Media Type	Thermal transfer ribbon, 1in wide roll
Media Capacity	1 roll, approximately 12,000 slides
Media Color	Black
Throughput	5 Slides/minute typical

SlideMate Available Accessories

Description	Order Number
Digital Imager	B35310006
19" All-in-one PC	A79510021
LabWriter Software for SlideMate	B81310001
Network Cable	B35310012
Crossover Ethernet Cable	B81310005
Keyboard Interface	B81320043
Keyboard	B81310004
Mouse	B81310019
Thermal Transfer Ribbon (hot-foil tape)	B81320041
Mini PC Package with built-in digital imager and clamp-style arm	B81300028
Mini PC Package with built-in digital imager and desktop stand	B81300029
Mini PC Package with clamp-style arm	B81300030
Mini PC Package with desktop stand	B81300031
Stylus, 3-pack	B81310022

To learn more about on-demand slide printing
Call: 800-522-7270 or
Visit: www.thermoscientific.com/slidemate

© 2010 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. ThinPrep is a registered trademark of Hologic Inc. All rights reserved.

Anatomical Pathology 4481 Campus Drive Kalamazoo, MI 49008 +1 (800) 522-7270

M41005 6/10

www.thermoscientific.com/pathology

Thermo
SCIENTIFIC